
Product Name: HSV Tag Mouse Monoclonal Antibody**Catalog #: AMM80045**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects HSV-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	HSV
Immunogen	Synthetic Peptide of HSV-Tag

Background

The HSV (herpes simplex virus) epitope tag is frequently engineered onto the N- or C- terminus of a protein of interest so that the tagged protein can be analyzed and visualized using immunochemical methods. The recognized HSV peptide epitope represents the amino acid sequence QPELAPEDPED.

Image Data



Western blot analysis of HSVtag in Multitag antigen 120ng, Multitag antigen 210ng lysates using HSVtag (6G5) antibody.