
Product Name: TAP Tag Mouse Monoclonal Antibody**Catalog #: AMM80047**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects TAP-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	TAP
Immunogen	Recombinant Protein of TAP Tag

Background

The TAP (Tandem Affinity Purification) method is an affinity purification method for the isolation of TAP-tagged proteins along with associated proteins. The TAP tag historically consists of a calmodulin binding peptide (CPB), a tobacco etch virus (TEV) protease cleavage site, and Protein A. However, additional tag combinations have been used with the TAP method including the combination of FLAG tags and HA tags.

Image Data



Western blot analysis of TAPtag (1F5) in TAP antigen 60ng, TAP antigen 100ng lysates using TAPtag (1F5) antibody.