

Product Name: SRT Tag Mouse Monoclonal Antibody**Catalog #: AMM80048**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects SRT-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	SRT
Immunogen	Recombinant Protein of SRT-Tag

Background

The SRT-tag is a peptide epitope comprised of 10 amino acids derived from the crystalline surface layer protein of *Rickettsia typhi*.

Image Data



Western blot analysis of SRTtag (5F5) in SRT antigen 60ng, SRT antigen 100ng lysates using SRTtag (5F5) antibody.