
Product Name: EYFP Tag Mouse Monoclonal Antibody**Catalog #: AMM80049**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects EYFP-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	EYFP
Immunogen	Recombinant Protein of EYFP

Background

Enhanced Yellow florescent protein. Yellow Fluorescent Protein (YFP) is a genetic mutant of green fluorescent protein, derived from *Aequorea victoria*. Its excitation peak is 514nm and its emission peak is 527nm.

Image Data



Western blot analysis of EYFPtag (3B2) in Multitag antigen 120ng, Multitag antigen 210ng lysates using EYFPtag (3B2) antibody.