

Product Name: ECFP Tag Mouse Monoclonal Antibody**Catalog #: AMM80052**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects ECFP-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------

Antigen Information

Gene Name	ECFP
Immunogen	Recombinant Protein of ECFP-Tag

Background

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Image Data

Western blot analysis of Recombinant ECFP Protein using ECFPTag Antibody.

