

Product Name: RFP Tag Mouse Monoclonal Antibody**Catalog #: AMM80055**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects RFP-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------


Antigen Information

Gene Name	RFP
Immunogen	Recombinant Protein of RFP-Tag

Background

RFP (Red Fluorescent Protein), cloned from *Discosoma coral* (DsRed, DsRed2 or drFP583), is a protein which exhibits bright red fluorescence when exposed to blue light.

Image Data



RFP fusion protein+ Primary antibody dilution at 1) 1:3000. 2) 1:5000 3) 1:10000