

Product Name: EBFP Tag Mouse Monoclonal Antibody**Catalog #: AMM80058**

For research use only.

Summary

Host	Mouse
Application	WB
Specificity	The antibody detects EBFP-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000
-----------------------	-------------------

Antigen Information

Gene Name	EBFP
Immunogen	Recombinant Protein of EBFP

Background

Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.