
Product Name: FITC-beta Tubulin Mouse Monoclonal Antibody**Catalog #: AMM80064**

For research use only.

Summary

Host	Mouse
Application	WB,ICC/IF
Specificity	endogenous
Reactivity	Human,Mouse,Rat,Goat,Hamster,Monkey
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:500-1:1000; IF-1:100-1:200
-----------------------	---------------------------------


Antigen Information

Gene Name	TUBB
Immunogen	Purified recombinant human Tubulin beta protein fragments expressed in E.coli.

Background

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Image Data



Immunocytochemistry analysis of beta Tubulin in HeLa cells using beta Tubulin antibody(FITC).