
Product Name: HRP-HA Tag Mouse Monoclonal Antibody**Catalog #: AMM80067**

For research use only.

Summary

Host	Mouse
Application	WB,ELISA
Specificity	The antibody detects HA-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG1
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS, containing 50% glycerol, pH7.4
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000; ELISA-1:5000-1:20000
-----------------------	---


Antigen Information

Gene Name	HA
Immunogen	Synthesized peptide of HA tag

Background

The HA tag is derived from the human influenza hemagglutinin corresponding to amino acids 98-106. It has been extensively used as a general epitope tag in expression vectors. This tag facilitates the detection, isolation, and purification of the proteins.

Image Data



Western blot analysis of HA Tag (3C6) in CHO-K1 (B) or CHO-K1 transfected with HAFusion proteins (A), using HATag antibody (HRP Conjugate).