

Product Name: HRP-GST Tag Mouse Monoclonal Antibody**Catalog #:** AMM80068

For research use only.

Summary

Host	Mouse
Application	WB,ELISA
Specificity	The antibody detects GST-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG2b
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS, containing 50% glycerol, pH7.4
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:10000; ELISA-1:5000-1:20000
-----------------------	---


Antigen Information

Gene Name	GST
Immunogen	Purified recombinant GST protein fragments expressed in E.coli.

Background

This antibody is high affinity to the glutathione-S-transferase (GST) from *Schistosoma japonicum*. This antibody will detect the expressed GST protein or a GST-fusion protein, and can be utilized for immuno-purification of GST-fusion proteins.

Image Data



Western blot analysis of purified recombinant GST fusion proteins with different molecular weights using GSTtag antibody (HRP Conjugate).