

Product Name: Flag-tag (M11) Rabbit Monoclonal Antibody**Catalog #: AMRe80001**

For research use only.

Summary


Host	Rabbit
Application	WB ,ELISA ,IF
Specificity	The antibody detects C-terminal, internal, and N-terminal Flag-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS,pH 7.4,containing 0.5%BSA and 50% Glycerol.
Purification	Protein A purified

Application


Dilution Ratio	WB 1:1000-5000, ELISA 1:20000-50000, ICC/IF 1:100-300
-----------------------	---

Antigen Information

Gene Name	Flag tag, Flag-tag,DDDDK TAG, DDDDK-TAG, DYKDDDDK tag,DYKDDDDK-tag
Immunogen	Synthetic Peptide of Flag-Tag

Background**Image Data**

Western Blot analysis of 293F and CHO cells transfected or non-transfected DDDK-tag expression vector by primary antibody at 1:3000 dilution.



Western Blot analysis of flag-tag protein using primary antibody at various dilution. Secondary antibody was diluted at 1:10000