

Production Name:Flag-tag (M11) Rabbit Monoclonal Antibody
Catalog#: AMRe80001



Summary

Production Name	Flag-tag (M11) Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB ,ELISA ,IF
specificity	The antibody detects C-terminal, internal, and N-terminal Flag-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG1
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS, pH 7.4, containing 0.5%BSAand 50% Glycerol.
Purification	Protein A purified

Immunogen Information

Gene Name	Flag tag; Flag-tag,DDDDK TAG, DDDDK-TAG, DYKDDDDK tag,DYKDDDDK-tag
Immunogen	Synthetic Peptide of Flag-Tag


Application

Dilution Ratio	WB 1:1000-5000, ELISA 1:20000-50000, IF-ICC 1:100-300
-----------------------	---

Background

Image Data

Production Name:Flag-tag (M11) Rabbit Monoclonal Antibody
Catalog#: AMRe80001



Western Blot analysis of 293F and CHO cells transfected or non-transfected DDDK-tag expression vector by primary antibody at 1:3000 dilution.

Note

For research use only .