

Product Name: HA tag Rabbit Monoclonal Antibody**Catalog #: AMRe80035**

For research use only.

Summary

Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP
Specificity	The antibody detects HA-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:500-1:1000; IHC-P-1:1000-1:5000; ICC/IF-1:500-1:1000; FC-1:200-1:500; IP-1:20-1:50
-----------------------	---

Antigen Information

Gene Name	HA
Immunogen	Synthesized peptide

Background

The HA tag is derived from the human influenza hemagglutinin corresponding to amino acids 98-106. It has been extensively used as a general epitope tag in expression vectors. This tag facilitates the detection, isolation, and purification of the proteins.