

Product Name: GST Tag Rabbit Monoclonal Antibody**Catalog #: AMRe80036**

For research use only.

Summary

Host	Rabbit
Application	WB,ICC/IF,IP
Specificity	The antibody detects GST-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:5000; ICC/IF-1:100; IP-1:20
-----------------------	---

Antigen Information

Gene Name	GST
Immunogen	Purified recombinant GST protein fragments expressed in E.coli.

Background

This antibody is high affinity to the glutathione-S-transferase (GST) from *Schistosoma japonicum*. This antibody will detect the expressed GST protein or a GST-fusion protein, and can be utilized for immuno-purification of GST-fusion proteins.