

Product Name: S Tag Rabbit Monoclonal Antibody**Catalog #: AMRe80044**

For research use only.

Summary

Host	Rabbit
Application	WB
Specificity	The antibody detects S-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:5000
-----------------------	------------------


Antigen Information

Gene Name	S
Immunogen	Synthesized peptide of S-tag

Background

S tag is an epitope tag composed of a 15 residue peptide, KETAAAKFERQHMDs, derived from the pancreatic ribonuclease A.

Image Data



Western blot analysis of S tag in multitag antigen 120ng, multitag antigen 210ng lysates using S Tag antibody.