

Product Name: MBP Tag Rabbit Monoclonal Antibody**Catalog #: AMRe80054**

For research use only.

Summary

Host	Rabbit
Application	WB
Specificity	endogenous
Reactivity	Escherichia Coli
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:5000
-----------------------	------------------


Antigen Information

Gene Name	malE
Immunogen	Recombinant Protein of MBP-Tag

Background

Part of the ABC transporter complex MalEFGK involved in maltose/maltodextrin import. Binds maltose and higher maltodextrins such as maltotriose.

Image Data



Western blot analysis of MBP lysates in E.coli using MBP antibody.