
Product Name: V5 Tag Rabbit Monoclonal Antibody**Catalog #: AMRe80060**

For research use only.

Summary

Host	Rabbit
Application	WB
Specificity	The antibody detects V5-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:1000-1:5000
-----------------------	------------------


Antigen Information

Gene Name	V5
Immunogen	Synthesized peptide of V5-tag

Background

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag antibody can be helpful in detecting the recombinant proteins, some of which include transmembrane and secreted proteins fusion protein.

Image Data



Western blot analysis of V5 tag in V5-tag fusion protein using V5 tag antibody.