

Product Name: HA-Tag Polyclonal Antibody**Catalog #: APRab80009**

For research use only.

Summary


Host	Rabbit
Application	WB, Elisa
Specificity	The antibody detects HA-tag fusion proteins.
Reactivity	Species independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.

Application

Dilution Ratio	WB 1:5000, ELISA 1:5000-20000, Optimal dilutions should be determined by the end user
-----------------------	---

Antigen Information

Gene Name	HA-tag
Immunogen	

Background**Image Data**

2ug HA fusion protein+ Primary antibody dilution at 1) 1:5000,2) 1:20000.
Secondary antibody was diluted at 1:20000