

Product Name: α -tubulin Polyclonal Antibody

Catalog #: APRab80018

For research use only.

Summary

Host	Rabbit
Application	WB, Elisa
Specificity	endogenous
Reactivity	Human, Mouse, Rat
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody detects endogenous α -tubulin protein.

Application


Dilution Ratio	WB 1:2000
-----------------------	-----------

Antigen Information

Gene Name	α -tubulin
Immunogen	Recombinant Protein of Tubulin alpha-1A chain

Background

Image Data



Western blot analysis of 1) Mouse skeletal muscle tissue, 2) Mouse kidney tissue, 3) Mouse Heart Tissue, 4) Mouse brain tissue, 5) Mouse Liver tissue, diluted at 1:2000. Secondary antibody was diluted at 1:20000