

Production Name: α -tubulin Polyclonal Antibody
Catalog#: APRab80018



Summary

Production Name	α -tubulin Polyclonal Antibody
Host	Rabbit
Application	WB, Elisa
specificity	
Reactivity	Human, Mouse, Rat
Isotype	
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 2% BSA and New type preservative N.
Purification	The antibody detects endogenous α -tubulin protein.

Immunogen Information


Gene Name	
Immunogen	Recombinant Protein of Tubulin alpha-1A chain

Application

Dilution Ratio	WB 1:2000
-----------------------	-----------

Background

Image Data



Western blot analysis of 1) Mouse skeletal muscle tissue, 2) Mouse kidney tissue, 3) Mouse Heart Tissue, 4) Mouse brain tissue, 5) Mouse Liver tissue, diluted at 1:2000.. Secondary antibody was diluted at 1:20000

Production Name: α -tubulin Polyclonal Antibody
Catalog#: APRab80018



Note

For research use only .