

Product Name: HRP-His Tag Rabbit polyclonal Antibody**Catalog #: APRab80070**

For research use only.

Summary

Host	Rabbit
Application	WB
Specificity	The antibody detects His-tag fusion proteins.
Reactivity	Species-independent
Isotype	IgG
Form	Liquid
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB-1:5000-1:10000
-----------------------	-------------------

Antigen Information

Gene Name	His
Immunogen	Synthesized peptide of His tag

Background

A polyhistidine-tag is an amino acid motif in proteins that consists of at least five histidine (His) residues, often at the N- or C-terminus of the protein. It is also known as hexa histidine-tag, 6xHis-tag, and by the trademarked name His-tag. Polyhistidine-tags are often used for affinity purification of polyhistidine-tagged recombinant proteins expressed in Escherichia coli and other prokaryotic expression systems.