

Summary

Production Name	FGFR4 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Ascitic Fluid

Immunogen

Gene Name	FGFR4
Alternative Names	FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD334
Gene ID	2264
SwissProt ID	P22455.

Application

Dilution Ratio	IF: 1:50-1:200
Molecular Weight	-

Background

Receptor for acidic fibroblast growth factor. Does not bind to basic fibroblast growth factor. Binds FGF19.


Product Name: FGFR4 Mouse Monoclonal Antibody
Catalog #: AMM00718




Research Area

Cardiovascular

Image Data



Western blot analysis of FGFR4 in Human FGFR4 (aa22369) lysates using FGFR4 antibody.



Immunofluorescence analysis of FGFR4 in HEK293 cells transfected with FGFR4hlgGfc using FGFR4 antibody (green) Blue: DRAQ5 fluorescent DNA dye.

Note

For research use only.