

**Product Name: IKB beta (5H3) Mouse Monoclonal Antibody**  
**Catalog #: AMM00735**

---

## Summary

<b>Production Name</b>	IKB beta (5H3) Mouse Monoclonal Antibody
<b>Description</b>	Mouse Monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	IHC-P
<b>Reactivity</b>	Human,Rat,Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	NFKBIB NFKBIB; IKBB; TRIP9; NF-kappa-B inhibitor beta; NF-kappa-BIB; I-kappa-B-beta; Ikb-B;
<b>Alternative Names</b>	Ikb-beta; IkappaBbeta; Thyroid receptor-interacting protein 9; TR-interacting protein 9; TRIP-9
<b>Gene ID</b>	4793
<b>SwissProt ID</b>	Q15653.

## Application

<b>Dilution Ratio</b>	IHC: 1:50-1:100
<b>Molecular Weight</b>	-

**Product Name: IKB beta (5H3) Mouse Monoclonal Antibody**  
**Catalog #: AMM00735**

---




## Background

IKB-beta is a protein of the NF-kappa-B inhibitor family. It inhibits NF-kappa-B by complexing with and trapping it in the cytoplasm. However, the unphosphorylated form resynthesized after cell stimulation is able to bind NF-kappa-B allowing its transport to the nucleus and protecting it from further IKBA-dependent inactivation.


## Research Area

Signal Transduction

## Image Data



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma using IKB beta (5H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IKB beta (5H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## Note

For research use only.