

Product Name: Prealbumin (1F6) Mouse Monoclonal Antibody
Catalog #: AMM00893



Summary

Production Name	Prealbumin (1F6) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	TTR
Alternative Names	Prealbumin; TBPA; Amyloidosis I; Transthyretin
Gene ID	7276
SwissProt ID	P02766.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background

Product Name: Prealbumin (1F6) Mouse Monoclonal Antibody
Catalog #: AMM00893




Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Research Area

Neuroscience


Image Data



Immunohistochemistry analysis of paraffin-embedded Human Liver Carcinoma Tissue using Prealbumin (1F6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Prealbumin (1F6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Product Name: Prealbumin (1F6) Mouse Monoclonal Antibody
Catalog #: AMM00893



Western blot analysis of Prealbumin (1F6) in Hela, Jurkat, mouse brain, MOLT4, Human serum lysates using Prealbumin (1F6) antibody

Note

For research use only.