

Product Name: FABP3 (3E8B10) Mouse Monoclonal Antibody
Catalog #: AMM03312



Summary

Production Name	FABP3 (3E8B10) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	FABP3 422 protein; Cardiac Fatty Acid Binding Protein; FABP 11; FABP 3; FABP11; FABP3; FABPH_HUMAN; fatty acid binding protein 11; Fatty acid binding protein 3; Fatty acid binding protein 3 muscle and heart; Fatty acid binding protein 3 muscle and heart mammary derived growth inhibitor; Fatty acid binding protein 3 muscle; Fatty acid binding protein 3
Alternative Names	
Gene ID	2170
SwissProt ID	P05413.

Application

Dilution Ratio	ELISA: 1:10000
Molecular Weight	-

Product Name: FABP3 (3E8B10) Mouse Monoclonal Antibody
Catalog #: AMM03312




Background

FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.

Research Area

Cardiovascular

Image Data



Standard Curve for HFABP: Capture Antibody mouse mAb (5C4F5) to HFABP at 4ug/ml and Detector Antibody mouse mAb(3E8B10u30017B10G11u30014H1B8)to HFABP at 0.08 ug/ml.

Note

For research use only.