

Product Name: TRIM33 (3H7) Mouse Monoclonal Antibody
Catalog #: AMM03434

Summary

Production Name	TRIM33 (3H7) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	TRIM33
Alternative Names	ECTO; PTC7; RFG7; TF1G; TIF1G; TIFGAMMA; TIF1GAMMA
Gene ID	51592
SwissProt ID	Q9UPN9.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 123 kDa; Observed MW: 140 kDa

Background

Acts as an E3 ubiquitin-protein ligase.


Product Name: TRIM33 (3H7) Mouse Monoclonal Antibody
Catalog #: AMM03434



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TIF1 gamma / TRIM33 in MCF-7 and HeLa lysates using TIF1 gamma / TRIM33 antibody.

Note

For research use only.