

**Product Name: PP1C gamma (2F2A3) Mouse
Monoclonal Antibody
Catalog #: AMM03442**

Summary

Production Name	PP1C gamma (2F2A3) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	PPP1CC PP 1G; PP-1G; PP1G; PP1G_HUMAN; PP1gamma; PPP 1G; PPP1CC; PPP1CC protein; PPP1G; Protein phosphatase 1 catalytic subunit gamma isoform; Protein phosphatase 1C catalytic subunit; Protein phosphatase 1C subunit; Protein phosphatase 2C gamma isoform; Serine/threonine phosphatase 1 gamma; Serine/threonine protein phosphatase PP1 gamma catalytic subunit; Serine/threonine-protein phosphatase PP1- gamma catalytic subunit.
Alternative Names	
Gene ID	5501
SwissProt ID	P36873.

Application

Dilution Ratio	WB: 1:500-1:1000
-----------------------	------------------

**Product Name: PP1C gamma (2F2A3) Mouse
Monoclonal Antibody
Catalog #: AMM03442**

Molecular Weight

Calculated MW: 37 kDa; Observed MW: 37 kDa


Background

The protein encoded by this gene belongs to the protein phosphatase family, PP1 subfamily. PP1 is an ubiquitous serine/threonine phosphatase that regulates many cellular processes, including cell division. It is expressed in mammalian cells as three closely related isoforms, alpha, beta/delta and gamma, which have distinct localization patterns. This gene encodes the gamma isozyme. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Research Area

Cell Biology

Image Data



Western blot analysis of PP1C in Hela, 293T and Jurkat lysates using PP1C antibody.

Note

For research use only.