

Product Name: SMARCC1 (2A4) Mouse Monoclonal Antibody
Catalog #: AMM03479



Summary

Production Name	SMARCC1 (2A4) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	SMARCC1
Alternative Names	BAF155; CRACC1; Rsc8; SMARCC1; SRG3; SWI3
Gene ID	6599
SwissProt ID	Q92922.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 123 kDa; Observed MW: 155 kDa

Background

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-

Product Name: SMARCC1 (2A4) Mouse Monoclonal Antibody
Catalog #: AMM03479




nucleosome topology). May stimulate the ATPase activity of the catalytic subunit of the complex.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of SMARCC1 in K562, Jurkat and HeLa lysates using SMARCC1 antibody.

Note

For research use only.