

**Product Name: CD44 (2H12) Mouse Monoclonal Antibody**  
**Catalog #: AMM03497**

---

## Summary

<b>Production Name</b>	CD44 (2H12) Mouse Monoclonal Antibody
<b>Description</b>	Mouse Monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	WB,IHC-F,IHC-P,ICC/IF,IP
<b>Reactivity</b>	Human

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	CD44 CD44; LHR; MDU2; MDU3; MIC4; CD44 antigen; CDw44; Epican; Extracellular matrix
<b>Alternative Names</b>	receptor III; ECMR-III; GP90 lymphocyte homing/adhesion receptor; HUTCH-I; Heparan sulfate proteoglycan; Hermes antigen; Hyaluronate receptor; Phagocytic glycopr
<b>Gene ID</b>	960
<b>SwissProt ID</b>	P16070.

## Application

<b>Dilution Ratio</b>	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
<b>Molecular Weight</b>	Calculated MW: 82 kDa; Observed MW: 82 kDa

**Product Name: CD44 (2H12) Mouse Monoclonal Antibody**  
**Catalog #: AMM03497**




## Background

Receptor for hyaluronic acid (HA). Mediates cell-cell and cell-matrix interactions through its affinity for HA, and possibly also through its affinity for other ligands such as osteopontin, collagens, and matrix metalloproteinases (MMPs). Adhesion with HA plays an important role in cell migration, tumor growth and progression.


## Research Area

Immunology

## Image Data




Western blot analysis of CD44 in HUVEC lysates using CD44 antibody.




Immunoprecipitation analysis of CD44 (2H12) in HeLa lysates using CD44 antibody.

**Product Name: CD44 (2H12) Mouse Monoclonal Antibody**  
**Catalog #: AMM03497**

---



Immunofluorescence analysis of CD44 (2H12) in HeLa cells using CD44 antibody(green),and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD44 (2H12) antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## **Note**

For research use only.