

Product Name: SMYD4 (7C8) Mouse Monoclonal Antibody
Catalog #: AMM03508

Summary

Production Name	SMYD4 (7C8) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	SMYD4
Alternative Names	ZMYND21.
Gene ID	114826
SwissProt ID	Q8IYR2.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 89 kDa; Observed MW: 89 kDa

Background


Product Name: SMYD4 (7C8) Mouse Monoclonal Antibody
Catalog #: AMM03508

Plays a critical role in cardiac development.


Research Area

Cancer

Image Data



Western blot analysis of SMYD4 in HeLa, MCF-7 and A431 lysates using SMYD4 antibody.



Western blot analysis of SMYD4 (7C8) in HeLa, MCF-7 and MDA-MB-468 lysates using SMYD4 antibody.

Note

For research use only.