

Product Name: DDX3 (6G8) Mouse Monoclonal Antibody
Catalog #: AMM03525

Summary

Production Name	DDX3 (6G8) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB, ICC/IF, IP, ChIP
Reactivity	Human, Rat, Mouse, Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	DDX3X
Alternative Names	DBX; DDX3; HLP2; DDX14.
Gene ID	1654
SwissProt ID	O00571.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200 IP: 1:20 ChIP: 1:20
Molecular Weight	Calculated MW: 73 kDa; Observed MW: 73 kDa

Background

Product Name: DDX3 (6G8) Mouse Monoclonal Antibody
Catalog #: AMM03525




Multifunctional ATP-dependent RNA helicase.


Research Area

Epigenetics and Nuclear Signaling


Image Data



Western blot analysis of DDX3 in HeLa, 3T3, C6, COS7, K562 and Jurkat lysates using DDX3 antibody.



Immunocytochemistry analysis of DDX3 (6G8) in HeLa using DDX3 antibody.



Immunoprecipitation analysis of DDX3 (6G8) in HeLa lysates using DDX3 antibody.



Product Name: DDX3 (6G8) Mouse Monoclonal Antibody
Catalog #: AMM03525

Note

For research use only.