

Product Name: ntpII (4E9) Mouse Monoclonal Antibody
Catalog #: AMM03534



Summary

Production Name	ntpII (4E9) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	neo
Alternative Names	neomycin phosphotransferase II
Gene ID	61686716.0
SwissProt ID	P00552.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 29 kDa; Observed MW: 29 kDa

Background

Resistance to kanamycin, neomycin, paromomycin, ribostamycin, butirosin and gentamicin B.


Product Name: ntpII (4E9) Mouse Monoclonal Antibody
Catalog #: AMM03534



Research Area

Signal Transduction

Image Data



Western blot analysis of ntpII in nontransgenic Brassica napus L. and transgenic Brassica napus L. lysates using ntpII antibody.

Note

For research use only.