

Product Name: GAD65 (6D5) Mouse Monoclonal Antibody
Catalog #: AMM03558

Summary

Production Name	GAD65 (6D5) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	GAD2
Alternative Names	GAD65
Gene ID	2572
SwissProt ID	Q05329.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 65 kDa; Observed MW: 60 kDa

Background

Catalyzes the production of GABA.


Product Name: GAD65 (6D5) Mouse Monoclonal Antibody
Catalog #: AMM03558



Research Area

Neuroscience

Image Data



Western blot analysis of GAD2 in rat Hippocampus and rat Brain lysates using GAD2 antibody.

Note

For research use only.