

**Product Name: DiMethyl-Histone H3 (Lys9) (10B2)**  
**Mouse Monoclonal Antibody**  
**Catalog #: AMM03582**

---

## Summary

<b>Production Name</b>	DiMethyl-Histone H3 (Lys9) (10B2) Mouse Monoclonal Antibody
<b>Description</b>	Mouse Monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	WB
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Methylated
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	H3C1
<b>Alternative Names</b>	H3K9me2; H3 histone; HIST1H3A; Histone cluster 1; H3a
<b>Gene ID</b>	8350
<b>SwissProt ID</b>	P68431.

## Application

<b>Dilution Ratio</b>	WB: 1:500-1:1000
<b>Molecular Weight</b>	Calculated MW: 15 kDa; Observed MW: 15 kDa

## Background

---

**Product Name: DiMethyl-Histone H3 (Lys9) (10B2)**  
**Mouse Monoclonal Antibody**  
**Catalog #: AMM03582**




H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.


## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of DiMethyl-Histone H3 in HeLa, rat Heart , Raw264.7 lysates using DiMethyl-Histone H3 (Lys9) (10B2) antibody



Western blot analysis of DiMethyl-Histone H3 (Lys9) (10B2) in MCF-7, C2C12, C6, 3T3, HeLa lysates using DiMethyl-Histone H3 (Lys9) (10B2) antibody.

## Note

For research use only.