

**Product Name: Flotillin 1 (6H9) Mouse Monoclonal Antibody**  
**Catalog #: AMM03590**

---



## Summary

<b>Production Name</b>	Flotillin 1 (6H9) Mouse Monoclonal Antibody
<b>Description</b>	Mouse Monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	WB
<b>Reactivity</b>	Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	FLOT1
<b>Alternative Names</b>	Flotillin-1
<b>Gene ID</b>	10211
<b>SwissProt ID</b>	O75955.

## Application

<b>Dilution Ratio</b>	WB: 1:500-1:1000
<b>Molecular Weight</b>	Calculated MW: 47 kDa; Observed MW: 47 kDa

## Background

---

**Product Name: Flotillin 1 (6H9) Mouse Monoclonal Antibody**  
**Catalog #: AMM03590**




May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.


## Research Area

Cell Biology

## Image Data



Western blot analysis of Flotillin1 in rat Brain, mouse Brain lysates using Flotillin1 (6H9) antibody.



Western blot analysis of Flotillin 1 (6H9) in mouse brain, rat brain lysates using Flotillin1 (6H9) antibody.

## Note

For research use only.