

Product Name: Transcription Initiation Factor IIB (2F6)
Mouse Monoclonal Antibody
Catalog #: AMM03639



Summary

Production Name	Transcription Initiation Factor IIB (2F6) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human,Mouse,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	GTF2B
Alternative Names	General transcription factor IIB; General transcription factor TFIIB; gtf2b; RNA polymerase II transcription factor IIB; S300 II; S300-II; TF IIB; TF2B; TF2B_HUMAN; TFIIB; Transcription initiation factor IIB.
Gene ID	2959
SwissProt ID	Q00403.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

Product Name: Transcription Initiation Factor IIB (2F6)
Mouse Monoclonal Antibody
Catalog #: AMM03639




Background

General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TFIIIB in Jurkat, MCF-7, COS7, 3T3 and Hela lysates using TFIIIB antibody.

Note

For research use only.