

Product Name: GRB2 (5G3) Mouse Monoclonal Antibody
Catalog #: AMM03656

Summary

Production Name	GRB2 (5G3) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purification

Immunogen

Gene Name	GRB2
Alternative Names	GRB2; ASH; Growth factor receptor-bound protein 2; Adapter protein GRB2; Protein Ash; SH2/SH3 adapter GRB2
Gene ID	2885
SwissProt ID	P62993.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

Background


Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.

Product Name: GRB2 (5G3) Mouse Monoclonal Antibody
Catalog #: AMM03656

Research Area

Signal Transduction

Image Data



Western blot analysis of GRB2 in K562, Raji, Jurkat, Ramos, C6, 3T3 and HeLa lysates using GRB2 antibody.

Note

For research use only.