

Product Name: ZNF364 Rabbit Monoclonal antibody
Catalog #: AMRe01416



Summary

Production Name	ZNF364 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	RNF115
Alternative Names	E3 ubiquitin-protein ligase RNF115; Rabring 7; Zinc finger protein 364
Gene ID	27246
SwissProt ID	Q9Y4L5.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 34 kDa

Background

Product Name: ZNF364 Rabbit Monoclonal antibody
Catalog #: AMRe01416




E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates and may play a role in diverse biological processes.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ZNF364/RFN115 in PC-12 lysates using ZNF364 antibody.



Western blot analysis of ZNF364/RFN115 in Jurkat, Hela lysates using ZNF364/RFN115 antibody.

Note

For research use only.