

Product Name: TNFAIP8 Rabbit Monoclonal Antibody
Catalog #: AMRe01424



Summary

Production Name	TNFAIP8 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	TNFAIP8
Alternative Names	TNF alpha-induced protein 8; NDED; SCC-S2
Gene ID	25816
SwissProt ID	O95379.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 21-23 kDa

Background

Acts as a negative mediator of apoptosis and may play a role in tumor progression.


Product Name: TNFAIP8 Rabbit Monoclonal Antibody
Catalog #: AMRe01424



Research Area

Immunology

Image Data



Western blot analysis of TNFAIP8 in CEM, K562 lysates using TNFAIP8 antibody.

Note

For research use only.