

Product Name: BAG2 Rabbit Monoclonal antibody
Catalog #: AMRe01434



Summary

Production Name	BAG2 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	BAG2
Alternative Names	BAG 2; Bag2
Gene ID	9532
SwissProt ID	O95816.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 24 kDa; Observed MW: 24 kDa

Background

Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release.


Product Name: BAG2 Rabbit Monoclonal antibody
Catalog #: AMRe01434



Research Area

Cell Biology

Image Data



Western blot analysis of BAG2 in Jurkat, HeLa lysates using BAG2 antibody.

Note

For research use only.