

Product Name: ARFRP1 Rabbit Monoclonal Antibody
Catalog #: AMRe01438



Summary

Production Name	ARFRP1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	ARFRP1
Alternative Names	ARP; Arp1; ARL18
Gene ID	10139
SwissProt ID	Q13795.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Trans-Golgi-associated GTPase that regulates protein sorting.


Product Name: ARFRP1 Rabbit Monoclonal Antibody
Catalog #: AMRe01438



Research Area

Signal Transduction

Image Data



Western blot analysis of ARFRP1/ARP in C6 lysates using ARFRP1 antibody.

Note

For research use only.