

Product Name: Desmocollin 3 Rabbit Monoclonal Antibody
Catalog #: AMRe01461



Summary

Production Name	Desmocollin 3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	DSC3
Alternative Names	Cadherin family member 3; Desmocollin-4; HT-CP
Gene ID	1825
SwissProt ID	Q14574.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 100 kDa; Observed MW: 100 kDa

Product Name: Desmocollin 3 Rabbit Monoclonal Antibody
Catalog #: AMRe01461


Background

Component of intercellular desmosome junctions.


Research Area

Signal Transduction

Image Data



Immunocytochemistry analysis of Desmocollin 3 (green) in HeLa using Desmocollin 3 antibody, and DAPI (blue).



Western blot analysis of Desmocollin 3 in K562, HeLa lysates using Desmocollin 3 antibody.

Note

For research use only.