

Product Name: PREB Rabbit Monoclonal Antibody
Catalog #: AMRe01472



Summary

Production Name	PREB Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	PREB
Alternative Names	Prolactin regulatory element-binding protein; Mammalian guanine nucleotide exchange factor mSec12
Gene ID	10113
SwissProt ID	Q9HCU5.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 45 kDa

Product Name: PREB Rabbit Monoclonal Antibody
Catalog #: AMRe01472




Background

Guanine nucleotide exchange factor that specifically activates the small GTPase SAR1B.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of PREB in K562, C6, 3T3, HeLa lysates using PREB antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PREB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.