

Product Name: GAD67 Rabbit Monoclonal Antibody
Catalog #: AMRe01511



Summary

Production Name	GAD67 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	GAD1
Alternative Names	Glutamate decarboxylase 1; 67 kDa glutamic acid decarboxylase; GAD 67; GAD1; GAD; GAD67
Gene ID	2571
SwissProt ID	Q99259.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 67 kDa; Observed MW: 67 kDa

Background

Product Name: GAD67 Rabbit Monoclonal Antibody
Catalog #: AMRe01511




Catalyzes the production of GABA.

Research Area

Signal Transduction

Image Data



Western blot analysis of GAD67 in HeLa lysates using GAD67 antibody.

Note

For research use only.