

**Product Name: Phospho-S6K2 (Ser371) Rabbit
Monoclonal Antibody
Catalog #: AMRe01555**

Summary

Production Name	Phospho-S6K2 (Ser371) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	RPS6KB2 RPS6KB2; STK14B; Ribosomal protein S6 kinase beta-2; S6K-beta-2; S6K2; 70 kDa
Alternative Names	ribosomal protein S6 kinase 2; P70S6K2; p70-S6K 2; S6 kinase-related kinase; SRK; Serine/threonine-protein kinase 14B; p70 ribosomal S6 kinase beta; S6K-beta; p70
Gene ID	6199
SwissProt ID	Q9UBS0.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 60 kDa

Product Name: Phospho-S6K2 (Ser371) Rabbit Monoclonal Antibody
Catalog #: AMRe01555


Background

Phosphorylates specifically ribosomal protein S6.

Research Area

Signal Transduction

Image Data



Western blot analysis of Phospho-p70 S6 Kinase beta/SRK (Ser371) in K562, rat Brain, C6, 3T3 lysates using Phospho-S6K2 (Ser371) antibody.

Note

For research use only.