

Product Name: Angiotensinogen Rabbit Monoclonal Antibody
Catalog #: AMRe01648

Summary

Production Name	Angiotensinogen Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	Agt
Alternative Names	Angiotensinogen; Serpin A8; AGT; SERPINA8; Des-Asp[1]-angiotensin II; Serine (or cysteine) proteinase inhibitor
Gene ID	11606.0
SwissProt ID	P11859.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Product Name: Angiotensinogen Rabbit Monoclonal Antibody
Catalog #: AMRe01648




Background

Essential component of the renin-angiotensin system (RAS), a potent regulator of blood pressure, body fluid and electrolyte homeostasis.

Research Area

Cardiovascular

Image Data



Western blot analysis of Angiotensinogen in mouse liver lysates using Angiotensinogen antibody.

Note

For research use only.