

Product Name: Annexin VI Rabbit Monoclonal Antibody
Catalog #: AMRe01651



Summary

Production Name	Annexin VI Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	ANXA6
Alternative Names	p68; p70; ANX6; CBP68; CPB-II
Gene ID	309
SwissProt ID	P08133.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 76 kDa; Observed MW: 70 kDa

Background

Product Name: Annexin VI Rabbit Monoclonal Antibody
Catalog #: AMRe01651




May associate with CD21. May regulate the release of Ca²⁺ from intracellular stores. Miscellaneous Seems to bind one calcium ion with high affinity.


Research Area

Cell Biology

Image Data



Western blot analysis of Annexin VI in K562, rat Brain, C6, 3T3 lysates using Annexin VI antibody.



Immunocytochemistry analysis of Annexin VI (green) in HeLa using Annexin VI antibody, and DAPI (blue)

Note

For research use only.