

Product Name: AS160 Rabbit Monoclonal Antibody
Catalog #: AMRe01681



Summary

Production Name	AS160 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	TBC1D4
Alternative Names	AS160; NIDDM5
Gene ID	9882
SwissProt ID	O60343.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 147 kDa; Observed MW: 160 kDa

Background

May act as a GTPase-activating protein for RAB2A, RAB8A, RAB10 and RAB14. Isoform 2 promotes insulin-induced glucose

Product Name: AS160 Rabbit Monoclonal Antibody
Catalog #: AMRe01681




transporter SLC2A4/GLUT4 translocation at the plasma membrane, thus increasing glucose uptake.

Research Area

Signal Transduction

Image Data



Western blot analysis of AS160 in 293 lysates using AS160 antibody.

Note

For research use only.