

Product Name: C9orf167 Rabbit Monoclonal Antibody
Catalog #: AMRe01747



Summary

Production Name	C9orf167 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	TOR4A
Alternative Names	C9orf167
Gene ID	54863
SwissProt ID	Q9NXH8.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Background

Predicted to enable ATP binding activity. Predicted to be located in extracellular region and platelet alpha granule lumen.

Product Name: C9orf167 Rabbit Monoclonal Antibody
Catalog #: AMRe01747




Predicted to be active in endoplasmic reticulum lumen and nuclear envelope.

Research Area

Cell Biology

Image Data



Western blot analysis of C9orf167 in K562, HeLa lysates using C9orf167 antibody.

Note

For research use only.