

**Product Name: CD134/OX40L Receptor Rabbit
Monoclonal Antibody
Catalog #: AMRe01776**

Summary

Production Name	CD134/OX40L Receptor Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	Tnfrsf4
Alternative Names	OX40 antigen; OX40L receptor
Gene ID	22163.0
SwissProt ID	P47741.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 30 kDa; Observed MW: 45 kDa

Background


**Product Name: CD134/OX40L Receptor Rabbit
Monoclonal Antibody
Catalog #: AMRe01776**

Receptor for TNFSF4/OX40L/GP34. Is a costimulatory molecule implicated in long-term T-cell immunity (By similarity).

Research Area

Immunology

Image Data



Western blot analysis of CD134 in mouse spleen lysates using CD134/OX40L Receptor antibody.

Note

For research use only.